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**GB** **Instruction for use**

NOBAALGIN® dressings are sterile dressings which consist of soft and highly absorbent alginate fibres. The non-woven structure of the dressing combined with its physical and chemical properties result in a high absorbency – particularly important in the treatment of moderate to heavily exuding wounds. The material adapts remarkably to the contours of the wound preventing it from drying out.

As soon as NOBAALGIN® dressings get in touch with wound exudate, an ion exchange process takes place. The firm fibre is modified into a moist gel. This transformation assists in maintaining a moist wound environment that supports the wound healing process and aids autolytic debridement. The gel can be removed atraumatically by means of forceps or glove-covered fingers.

**Product composition:**  
Alginate Dressing composed of calcium alginate fiber and polysorbate.

**Intended use:**  
Alginate dressing is for all moderate to heavily exuding or bleeding wounds.

**Notes:**

- Before using a sterile product, visually inspect the packaging to ensure that it is intact. Do not use if individual package is damaged or open.

- Single use device: The reuse of single use devices may lead to a microbiological hazard. Treatment for reuse may adversely affect the performance of the product.

Should serious incidents occur in relation to the device, they are, according to MDR (EU) 2017/745, to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Precaution:**  
Medical assessment of the condition of the wound and the causes of the delayed wound healing are necessary before treating wounds with impaired healing. Treatment with the product cannot replace treatment of the cause of the impaired healing. In the event of full thickness second degree and third degree burns, product may be used under medical supervision following debridement of necrotic tissue or surgical treatment.

**Directions for use:**

- Cleanse the wound in accordance with standard nursing practice.
- Apply the dressing to the wound and secure it with a suitable secondary dressing.
- In case of dry wounds the product can be moistened with normal saline.
- The Product may be left in-situ for extended periods (up to several days), depending on the level of wound secretion. The dressing should be changed once it has fully turned into a gel in accordance with nursing standards.
- Should any fibres remain in the wound, these may be easily removed by a wound irrigation solution and are broken down naturally by the body.

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