



NOBAMED Paul Danz AG

RUDATAPE®

REF 070003

Product Description, Intended use, Application

The non-elastic white tape (3.8 cm x 10 m) is provided with a zinc oxide adhesive which develops its highest adhesive strength at body temperature. Owing to its special fabric structure, RUDATAPE® is easy to tear. It is used for functional tape bandages for support and relief in the field of surgery and orthopedics. In case of wounds, they are covered with a suitable wound dressing first.

Composition

Viscose, zinc oxide adhesive
(colophoniumfree)

Contraindications

Do not use in case of extensive haematomas, extensive skin injuries and a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

RUDATAPE® is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box
made of cellulose

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and transport

To be stored in a dry and dust-free environment, protected against sunlight.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.