



**NOBAMED Paul Danz AG**

## ABDOMINALVERBAND

**REF 371200**

### Product Description, Intended use, Application

The abdominal bandage, size 4: 95-105 cm, is permanently elastic and 25 cm wide. It is provided with a padding which is to lie on the abdomen. Owing to a hook and loop fastening the bandage has an adjustment tolerance of approximately 10 cm. It serves to relieve the wound suture after operations in the abdominal cavity. The abdominal pad distributes the pressure on the abdominal wall evenly. The bandage is, fasteners closed, washable at 30° C in the case of domestic use. Due to its material composition, the product may not be treated chemo-thermally and may therefore not be used on different patients. The abdominal bandage is not intended for permanent use as a dorsal bandage.

### Composition

Cotton, polyester, elastane, polyamide

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary: folding box  
made of cellulose

Secondary: box  
made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

### Storage and Transport

Dry and dustfree.

### Single use device

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.