



**NOBAMED Paul Danz AG**

## NOBA<sup>®</sup> WASH CAP

**REF 677000** (as of LOT 136653)

### Product Description, Intended use, Application

This complete hair care in a disposable cap without rinsing is designed for the professional hair care in particular of bedridden patients. The mild and non-irritating formula treats and cleans hair without water. The cap can be warmed in its packaging in a microwave oven at 500-700 watts for 20 seconds. Test temperature before use. Cover hair with the cap, massage thoroughly, remove cap, dry hair as usual.

### Composition

Nonwoven, Aqua, Disodium Cocoamphodiacetat, Dimethicon, Glycerine, Phenoxyethanol, Hamamelis Virginiana Extract, Aloe Barbadensis Leaf, DMDM Hydantoin, Ethylhexylglycerin, Benzalkonium Chloride, Tocopheryl Acetate, Citric Acid, perfume.

### Warnings, Notes

Do not dispose of in toilet.  
Keep out of reach of children.  
Avoid contact with the eyes. In case of contact with the eyes, rinse immediately with plenty of water. If irritation persists, seek medical advice.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE-bag  
Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.