

PRODUCT SPECIFICATION

REF 1911 CEDERROTH 4-IN-1 MINI



1 GENERAL PRODUCT DESCRIPTION

Description

Cederroth 4-in-1 mini is a sterile universal dressing, which consists of one elastic bandage and one compress with a non-stick wound contact layer which is stitched to the bandage. The dressing is wrapped in a pouch with printed pictorial user instructions.

The Cederroth 4-in-1 mini dressing is intended to be used in first aid situations and it is especially suitable for fingers and toes. The dressing has four functions and thereby replaces a variety of traditional first aid dressings and compresses:

- Pressure pad for heavy bleeding
- Protective dressing for superficial wounds
- Burn dressing with special wound protecting surface
- Temporary support dressing

The product is intended for single use.

Storage information

Dry storage in room temperature.

For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (w x l mm)</u>	
		<u>Bandage</u>	<u>Compress</u>
Mini blood stopper in sterile package, packed in outer plastic pouch with illustrations	1	3000 x 60	120 x 80

Material

Bandage:	Viscose
Compresses:	Absorbent layer: Viscose staple fibers
	Non-adherent outer layer: Porous net made of high density polyethylene
Sterile barrier:	Paper and transparent film made of polyethylene/polyamide.
Outer pouch:	Polyethylene with paper insert.

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610019112	7310801911003	7312021911009
Content:	1 dressing	30 CU	35 SKU/layer x 6 layers
Size (w x l x h mm):	40 x 25 x 65	155 x 210 x 105	800 x 1200 x 1080
Gross weight:	14 g	520 g	134 kg
Packing material:	Plastic 2 g, Paper 2 g	Cardboard 82 g	-

