



NOBAMED Paul Danz AG

NOBANETZ®

REF 302005

Product Description, Intended use, Application

NOBANETZ® is a disposable tubular net which is soft, elastic and breathable. It is 25 m long and suitable for the toe (size 0.5). A colour-coding system (green) facilitates the distinction of sizes. The tubular net can be cut to length and is easy to handle. Due to the strong fabric the net can be cut to length without fraying. NOBANETZ® is used for the quick and easy fixation of wound pads and dressings. It is also used as a cover over bandage dressings or plaster splints. The product is individually packed.

Composition

Polyamide, spandex

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging:	foil bag
Secondary packaging:	folding box made of cellulose
Tertiary Packaging	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.