



NOBAMED Paul Danz AG

WÖCHNERINNENVORLAGE

REF 001034

Product Description, Intended use, Application

The absorbent, non-sterile maternity pad serves for the absorption of the lochia after the delivery. There is an adhesive strip on the back site. The size of the maternity pad is 8 cm x 34 cm.

Composition

Pulp, polypropylene

Residual Risks, Undesirable Side Effects, Warnings, Notes

The product must be sterilized in accordance with a validated sterilization method before use on open wounds.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.