



NOBAMED Paul Danz AG

RIBOMEN®

Product Description and Purpose

Absorbent incontinence pads for men suitable for slight to moderate urinary leakage. The pad is fixed with adhesive strips and provides an excellent fit. RIBOMEN® is anatomically shaped and has a leakage protection on both sides. A special moisture dispensation layer protects against re-wetting.

Composition

Cellulose, polypropylene, polyethylene, elastomers, polyacrylate adhesive

Normative and Legal Requirements

RIBOMEN® is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The absorption capacity according to ISO 11948 is 600 g.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

