



**NOBAMED Paul Danz AG**

## **RIBOSLIP®**

**Size M**

**REF 006220** (as of LOT 134549)

### **Product Description, Intended use, Application**

Body-shaped incontinence diapers for adults, for severe incontinence. RIBOSLIP® is gentle to the skin and is fixed by means of two closure stripes on each side. A special moisture dispensation layer protects against re-wetting. RIBOSLIP® has indicator strips.

### **Composition**

Cellulose, superabsorber, polypropylene, polyethylene, elastomers

### **Contraindications**

The product should not be used in the case of a known allergy against the material.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The absorption capacity according to ISO 11948 is higher than 2300 g.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary: PE bag

Secondary: carton made of cellulose

### **Symbols used in labelling**



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### **Storage and Transport**

To be stored in a dry and dust-free environment

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.