



NOBAMED Paul Danz AG

RIBOSLIP®

Size L

REF 006230 (as of LOT 134550)

Product Description, Intended use, Application

Body-shaped incontinence diapers for adults, for severe incontinence. RIBOSLIP® is gentle to the skin and is fixed by means of two closure stripes on each side. A special moisture dispensation layer protects against re-wetting. RIBOSLIP® has indicator strips.

Composition

Cellulose, superabsorber, polypropylene, polyethylene, elastomers

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The absorption capacity according to ISO 11948 is higher than 2300 g.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary: PE bag

Secondary: carton made of cellulose

Symbols used in labelling



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.