



NOBAMED Paul Danz AG

RIBOCARE®

REF 008104

Product Description and Purpose

RIBOCARE® light bed sheet (60 x 60 cm, 4 layers) consists of a liquid-proof foil several cellulose layers. The filling is covered by a non-woven which is closed on all sides. The bed sheet is used for the care of incontinence patients. In addition to that, it also serves as protection for the bed when moist wound dressings are used or dressings are changed. The cellulose layers safely absorb liquids. The foil prevents the liquids from discharging and so protects the bedlinen from becoming soiled. RIBOSOFT® supports the nursing staff in the hospital as well as in the domestic home care of outpatients.

Composition

Fresh cellulose, polypropylene, polyethylene

Normative and Legal Requirements

RIBOSOFT® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1.

The absorption capacity according to ISO 11948 is 350 g.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears the following marking and symbols

