



NOBAMED Paul Danz AG

RIBOPANTS®

Size L

Product Description and Purpose

Closed incontinence pants for adults, suitable for moderate to severe incontinence. RIBOPANTS® are worn like normal underpants. The pants are equipped with elastic cuffs at the leg cut-outs and the waist. A special moisture dispensation layer protects against re-wetting.

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

Composition

Cellulose, polypropylene, polyethylene, elastomers

The product bears following labels and symbols:

Normative and Legal Requirements

RIBOPANTS® is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The absorption capacity according to ISO 11948 is 910 g.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.



Packaging

Primary packaging: PE bag

Date of information: April 01, 2020 [REV 4]

replaced: September 25, 2015 [REV 3]