



**NOBAMED Paul Danz AG**

## **RUDADETEC®**

**REF 040307**

### **Product Description and Purpose**

Blue plaster strips (2.5cm x 7.6cm) with worked-in aluminium strips which render the plasters detectable - specifically produced for the food industry. For smaller, primary healing wounds.

### **Composition**

Polyethylene, polyacrylate adhesive, aluminium

### **Normative and Legal Requirements**

RUDADETEC® is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment.

### **The product bears the following marking and symbols**



0123