



NOBAMED Paul Danz AG

RUDANASAL®

REF 069250

Product Description and Purpose

Special plaster made of brown, elastic non-woven fabric for the fixation of nasally inserted tubes, like stomach tube, feeding tube, oxygen tube, endotracheal tube. Size: adults. The plaster is provided with an excellently adherent, hypoallergenic polyacrylate adhesive. Owing to its hypoallergenic features, RUDANASAL® fixation plasters may also be used for patients prone to allergies.

Composition

Polyester, viscose, polyacrylate adhesive

Normative and Legal Requirements

RUDANASAL® is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

