



NOBAMED Paul Danz AG

RUDAVLIES®

REF 070130

Product Description and Purpose

RUDAVLIES® is a transversely elastic soft and white nonwoven fabric coated with hypoallergenic polyacrylate adhesive. 10 m are on a roll (width 30 cm). RUDAVLIES® can be cut to size individually and is used to fix primary wound dressings. In addition, RUDAVLIES® can be used for the fixation of e.g. catheters, tubes and measuring elements.

Contra-Indications

Known allergy against one of the components listed below.

Composition

Polyester, polyacrylate adhesive

Normative and Legal Requirements

RUDAVLIES® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, according to rule 1

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	folding box made of cellulose
Secondary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

