



NOBAMED Paul Danz AG

RUDAVLIES®-steril

REF 072586

Product Description and Purpose

The carrier material of the sterile wound dressing (size: 10 cm x 6 cm) consists of a white, elastic non-woven and is coated with a hypoallergenic polyacrylate adhesive. The wound pad does not adhere to the wound and is completely protected by the carrier material. The dressings is individually packed in sterile packaging.

Contraindications

Known allergies to one of the ingredients listed below.

Composition

Polyester, viscose, polyethylene, polyacrylate adhesive

Normative and Legal Requirements

RUDAVLIES®-steril is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper- packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected from sunlight

The product bears following labels and symbols:



STERILE EO



CE
0123