



NOBAMED Paul Danz AG

RUDAFILM®

REF 074110

Product Description and Purpose

Transparent polyurethane film with polyacrylate adhesive. For fixation of primary dressings. RUDAFILM® is waterproof, impermeable to bacteria, but permeable to vapor. 10 m x 10 cm are supplied as a roll and can be individually cut to size. It is not to be used directly on injured skin.

Contra-Indications

Known allergy to one of the ingredients listed below.

Composition

Polyurethane, polyacrylat adhesive

Normative and Legal Requirements

RUDAFILM® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, according rule I

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	folding box made of cellulose
Secondary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

