



NOBAMED Paul Danz AG

RUDAMED®-sensitive

REF 082506

Product Description, Intended use, Application

The white adhesive dressing (6 cm x 5 m) is made of a transversely elastic, soft non-woven fabric and is provided with a hypoallergenic polyacrylate adhesive. The dressing is therefore particularly suitable for sensitive skin. RUDAMED®-sensitive is used to cover smaller primary healing wounds. The non-adherent wound pad absorbs blood and exudate. RUDAMED®-sensitive is packed in a convenient dispenser box with opening.

Composition

Polyester, viscose, polyethylene and polyacrylate adhesive

Contra indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.