



**NOBAMED Paul Danz AG**

## **RUDAMED®-Junior**

**REF 084106**

### **Product Description, Intended use, Application**

The carrier material of the adhesive dressing (6 cm x 10 cm) is made of a water-resistant, flexible synthetic printed with color designs for children. RUDAMED®-junior is used to cover smaller, primary healing wounds. The non-adherent wound pad absorbs blood and exsudate. Owing to its polyacrylate adhesive coating the product is suitable for sensitive skin.

### **Composition**

Polyethylene, polypropylene, viscose and polyacrylate adhesive

### **Contra indications**

The product should not be used in the case of a known allergy against the material.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Legal Requirements**

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

|            |                                     |
|------------|-------------------------------------|
| Primary:   | PE bag                              |
| Secondary: | folding box<br>made of<br>cellulose |
| Tertiary:  | carton made of<br>cellulose         |

### **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### **Storage and Transport**

Dry and dustfree

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.