



NOBAMED Paul Danz AG

RUDADERM®- Fingerkuppenpflaster

REF 089090

Product Description and Purpose

Individually packaged sterile plaster with flexible, textile carrier material. Due to the special shape smaller, primarily healing wounds on fingertips can be covered quickly and easily.

Contraindications

Known allergy to one of the ingredients listed below.

Composition

Polyester, viscose, polypropylene, polyethylene, hypoallergenic adhesive

Normative and Legal Requirements

The fingertip plaster is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class Is according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper-paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears the following symbols and labels:

