



**NOBAMED Paul Danz AG**

## **RUDADERM®- Zehenpflaster**

**REF 089100**

### **Product Description and Purpose**

Individually packaged sterile plaster with flexible, textile carrier material. Due to the special shape smaller, primarily healing wounds on toes can be covered quickly and easily.

### **Contraindications**

Known allergy to one of the ingredients listed below.

### **Composition**

Polyester, viscose, polypropylene, polyethylene, hypoallergenic adhesive

### **Normative and Legal Requirements**

The toe plaster is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class Is according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: paper-paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment

**The product bears the following symbols and labels:**

