



NOBAMED Paul Danz AG

RUDADERM[®]-junior

REF 091105

Product Description, Intended use, Application

The plaster stripes consist of a carrier material (1.6 x 5.7cm) which is made of a flexible, air-permeable plastic with colorful animal motifs and is provided with a hypoallergenic adhesive. The non-sticking wound pad (1 x 1.7 cm) is protected all around by the carrier material. The individually sterile packed strips are used to cover primarily healing smaller wounds, injection sites and small punctures (for example vein punctures).

Composition

Polyethylene, viscose, polyester, polypropylene, hypoallergenic adhesive

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging:	Paper packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	Carton made of Cellulose

Symbols used in labelling

Explanations at www.nobamed.com



STERILE EO



CE
0123

Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree.

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.