



NOBAMED Paul Danz AG

NOBAVIT®

REF 145030

Product Description, Intended use, Application

The disposable product NOBAVIT® is a white, soft, and absorbent gauze and cotton tissue roll (30 cm x 3 m). The roll consists of cotton wool wrapped in gauze. NOBAVIT® is closed at the edge and is suited for creating padded dressings of all kinds. The product is individually packed.

Composition

Cotton

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

Cotton wool according to BP or USP absorbent cotton

Gauze according to Din EN 14079

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.