



Katheter-Set, latexfrei

REF 190600

Contents of this procedure pack:

1 x kidney-shaped tray, NOBAMED
Paul Danz AG,

2 x gauze swab 7,5 x 7,5 cm, 8-ply,
NOBAMED Paul Danz AG,

6 x gauze ball 3,0 cm, NOBAMED
Paul Danz AG,

1 x NOBAWRAP® 60 x 60 cm wrap-
ping drape, NOBAMED Paul Danz AG,

2 x NOBAGLOVE® -Nitril size M,
NOBAMED Paul Danz AG,

2 x forceps, blue, NOBAMED Paul Danz
AG,

1 x NOBAWRAP® 60 x 60 cm,
NOBAMED Paul Danz AG,

1 x plastic tray 125, NOBAMED Paul
Danz AG,



This product is for single and tempo-
rary use on one patient and must be applied
by specially trained medical personnel. Any
information for use that is available for indi-
vidual products must be observed!

Normative and Legal Requirements

Procedure Pack according to Art. 12 of the
Medical Device Directive 93/42 / EEC.

The compatibility of the customized medical
devices or other products has been as-
sessed with due regard for the indications of

the respective manufacturers, and opera-
tions have been performed in accordance
with the aforementioned indications.

The manufacturers' intended use is not
changed.

The systems and procedure packs are
packed in compliance with a validated pro-
cedure.

All activity is properly monitored and con-
trolled as part of our QM system, in accord-
ance with DIN EN ISO 13485:2016.

A validated sterilization process is used for
sterile procedure packs, taking into account
the manufacturers' instructions, so that
product functionality and product safety of
the components are maintained after a sys-
tem or procedure pack has been compiled
and sterilized.

**Sterilization and packaging are subject
to supervision of our Notified Body TÜV
Süd PS 0123 pursuant to Annex V.**

The product does not contain dangerous
toxic substances according to REACH. It
has CE marking and DIN EN ISO 15223-1
labels on all its packaging

Storage

To be stored in a dry and dust-free envi-
ronment, protect from sunlight.

Symbols used in labelling

Explanation at www.nobamed.com

