



NOBAMED Paul Danz AG

NOBACOSMED®

REF 201100

Product Description and Purpose

White, disposable tissues made of soft cellulose, medically bleached. They can be used for the pre-cleaning of medical devices, for example for removing of ultrasound gel from ultrasound probes. Beside that there are diverse applications for the tissues such as using them as a wipe for cleaning the skin.

Composition

Cellulose

Normative and Legal Requirements

NOBACOSMED® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, according to rule I

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary-
packaging:

Folding box made
of cellulose

Secondary-
packaging:

Carton made
of cellulose

Lagerung

To be stored in a dry and dust-free environment

