



NOBAMED Paul Danz AG

NOBATISSUE®-weiß

REF 206534

Product Description, Intended use, Application

Non-sterile, **highly bleached** dressing cellulose, several individual and creped layers placed on top of each other, **size 28 cm x 37 cm**. They are absorbent, white, soft, and odourless. The highly bleached NOBATISSUE®-weiß products are mainly used in the health and care sector as absorbent material and for protection or padding purposes. NOBATISSUE®-weiß is for single use only.

Composition

Highly bleached (oxygen bleach) dressing cellulose

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDR (EU) 2017/745*.

The cellulose layers comply with the requirements for dressing cellulose of the German Pharmacopoeia, edition 10.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: Carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.