



NOBAMED Paul Danz AG

NOBAZELLTUPF®-steril

REF 209405

Product Description and Purpose

The sterile cellulose swab roll is provided with perforations, which allow an easy taking of the swabs. Each swab has a size of 4 cm x 5 cm.

The swabs are suited for the preparation of simple punctures with low risks, for example vein punctures for blood sampling.

Contraindications

Sterilized cellulose swabs are not suitable for the preparation of punctures with a higher risk, e.g. punctures of joints or organs. In this case individually packed sterile swabs are to be used.

Not to be used in case of a known allergy against the component listed below.

Composition

Cellulose

Normative and Legal Requirements

NOBAZELLTUPF®-steril is a medical product according to the MPG and the directive 93/42/EEC and is classified as class Is according to rule 4.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: film-film packaging

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

