



NOBAMED Paul Danz AG

NOBAPANTS®

REF 311015

Product Description, Intended use, Application

The disposable product NOBAPANTS® (size **XXXL**) is a white fixation panty made of thin spunbonded non-woven fabric, reinforced at the crotch. It has elastic leg and hip cuffs that ensure a good fit. NOBAPANTS® is used to secure absorbent pads and underpads. It is also frequently used during medical examinations and treatments to reduce the feeling of shame.

Composition

Polypropylene, elastomer

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE-bag

Secondary packaging: Carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.