



NOBAMED Paul Danz AG

ZWIRNHANDSCHUHE

REF 323130

Product Description, Intended use, Application

The white disposable knitted gloves, **size 8.5**, are made of medically bleached cotton. They are suitable for use in patient care and dermatology. They are mainly used to prepare the disinfected surgical instruments in sterile areas as they protect the instruments against the natural acid and fat content of the hands. Cotton gloves may also be worn underneath other gloves.

Composition

Cotton

Notes

The product has been subjected to a validated microbial reduction process and is therefore pre-sterilized.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PP-bag
cardboard

Secondary packaging: carton made of
cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.