



## WUNDHAKENBEZUG

REF 324034

### Product Description and Purpose

Disposable single-layer white covers sewn of medically bleached knitted cotton fabric, 10 cm x 26 cm, packed per 20 pcs, with worked in X-ray thread and with bands. The bands are used to fix the covers at the wound hook. The covers are put on the wound hooks before sterilization, to protect the instruments. The covers may be left on the wound hooks during their use and will thus reduce light reflections caused by the metal due to the operation lamp. At the same time sensitive tissue will be protected against the edges of the wound hooks. Before use on open wounds, the product must be sterilized in accordance with a validated sterilization method (Moist Heat in accordance with DIN EN ISO 17665, Radiation in accordance with DIN EN ISO 11137, or Ethylene Oxide in accordance with DIN EN ISO 11135).

### Composition

Cotton,  
X-Ray thread contains barium sulphate,  
PVC and DINP

### Normative and Legal Requirements

WUNDHAKENBEZUG (Wound Hook Cover) is a medical product according to the MDD 93/42/EEC and the MDR (EU) 2017/745 and is classified as class IIa, rule 7.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

### Storage

To be stored in a dry environment, protected against direct sunlight

The product bears the following symbols and marking:

