



**NOBAMED Paul Danz AG**

## NOBA<sup>®</sup> HEMD

**REF 324044**

### Product Description, Intended use, Application

The disposable product NOBA<sup>®</sup> HEMD, **for adults, small (for women up to size 44 or men up to size 50)**, is a finely woven, very tightly knitted shirt, natural- coloured, with a ribbed structure. Its high elasticity guarantees best adaptability. The unbleached cotton is water-repellent. Due to its natural fibre, the material has hypoallergenic features and is very skin-tolerant. The shirt is used under rigid or half-rigid cast bandages. In the case of dermatological diseases, the NOBA<sup>®</sup> HEMD serves as cover for ointment dressings or similar forms of therapy.

### Composition

Cotton

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: Folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.