



NOBAMED Paul Danz AG

CLAVICULA-Bandage

REF 372035

Product Description and Purpose

The one-piece special bandage (bambi – children's size 92- 104) is suited for the treatment of clavicle fractures.

It has a terry cloth cover on the interior side and is padded in the arm-pit area and under the tensioning rings. This increases wearing comfort. If the bandage is stretched sufficiently, the ends are fixed with the sewn-on hook and loop fasteners.

The bandage is, fasteners closed, washable at 30° C in the case of domestic use. Due to its material composition, the product may not be treated chemo-thermally and may therefore not be used on different patients.

Composition

Cotton, Polyester, Polyamide, Polyurethane

Normative and Legal Requirements

The Clavicle-Bandage is a medical product according to directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

