



**NOBAMED Paul Danz AG**

# FIXIERGURT

**REF 374675**

## **Product Description, Intended use, Application**

The padded fixation belt (7.5 cm x 6 m) is made of foam which is covered with a cotton tube. It can be used to make traditional figure-of-eight bandages or arm slings.

## **Composition**

Cotton, polyurethane, polyamide

## **Contraindications**

The product should not be used in the case of a known allergy against the material.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Normative and Legal Requirements**

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

|            |                                     |
|------------|-------------------------------------|
| Primary:   | folding box<br>made of<br>cellulose |
| Secondary: | carton made of<br>cellulose         |

## **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

## **Storage and Transport**

Dry and dust free

## **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

## **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.