



NOBAMED Paul Danz AG

NOBA[®]-KNIESCHIENE

Product Description and Purpose

NOBA[®]-KNIESCHIENE is a coat-shaped knee splint with a 15° flexion. It is used to immobilize the joint in the case of traumatic or degenerative changes of the knee joint, in the rehabilitation phase after knee operations or treatments with P.O.P. bandages, stabilization after knee joint endoprosthesis implantation, ruptures of the crucial ligament, complete instability and in the case of minor or medium range instability of the ligaments. The flexion of the splint can be adapted individually. The splint consists of three parts and can therefore be adapted to the different leg circumferences. The knee splint is washable at max. 30° C, hook and loop fasteners closed, if used domestically. Due to its material composition, the product may not be treated chemo-thermally and may therefore not be used on different patients.

Composition

Cotton, polyamide, polyester

Normative and Legal Requirements

NOBA[®]-KNIESCHIENE (Knee Splint) is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

