



NOBAMED Paul Danz AG

NOBAFIX®

Individually packed in foil

REF 410420

Product Description, Intended use, Application

This fixation bandage with a width of 12 cm and a stretched length of 4 m is a rolled, lengthwise-elastic dressing retention bandage. Due to its weave, the bandage has a fine, smooth and tight structure. The elastic bandage is used to fix wound pads, especially over joints and body curves. As a damp bandage, NOBAFIX® is also used to shape and/or fix plaster splints and cast dressings. NOBAFIX® is for single use.

Composition

Polyester

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PP foil

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree, protect from sunlight.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.