



NOBAMED Paul Danz AG

IDEALBINDE

Individually packed in foil

REF 447520

Product Description, Intended use, Application

The textile elastic, natural coloured and short stretched cotton crepe bandage is **20 cm wide and 5 m long** (in stretched condition). It has fast edges. Owing to its material and its weave, it is hypoallergenic and universally applicable in all medical fields. The rolled bandage is for single use only.

Composition

cotton

Contraindications

The product should not be used in the case of a known allergy against the material.

Note

Please ensure after application, that the bandage must neither cause pressure sores, strangulation marks nor pains, numbness and circulatory disorders.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The bandage is manufactured in accordance with the requirements of DIN 61632.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PP foil

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.