



**NOBAMED Paul Danz AG**

## NOBA® OP-SOCKEN

REF 660100

### Product Description, Intended use, Application

The disposable product NOBA® OP-SOCKEN are soft, elastic, pre-rolled socks in one size, without heel seams and with elastic cuffs. The white socks are suitable for preoperative preparation of patients. Socks contribute to the prevention of perioperative hypothermia by protecting the patient from heat loss as a preventive measure. In addition, NOBA® OP-SOCKEN are universally applicable and counteract the spread of germs in sensitive areas.

### Composition

Polyester, elastane

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: foil bag

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.