



**NOBAMED Paul Danz AG**

## **NOBASUIT®-Thermojacke**

**REF 663530 (as of LOT 135289)**

### **Product Description, Intended use, Application**

The disposable NOBASUIT® thermal jacket in blue/ green, **size XXXL**, serve as scrub clothing and clean air clothing to reduce the spreading of infectious agents. The material is anti-static, breathable, low-particle, convenient and comfortable to wear. The thermal jacket is provided with knitted cuffs at the neck area and arm ends. It has two pockets and is equipped with white buttons in the front.

### **Composition**

Polypropylene (SMS), polyester

### **Contraindication**

The product should not be used in the case of a known allergy against the material.

### **Note**

The product has been subjected to a validated microbial reduction process and is therefore pre-sterilized.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to MDR (EU) 2017/745.

The product complies with the requirements of DIN EN 13795-2 for clean air suits.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

### **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

### **Storage and Transport**

Dry and dustfree.

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.