



**NOBAMED Paul Danz AG**

## **BESUCHERKITTEL (Visitor's Gown)**

**REF 665130**

### **Product Description, Intended use, Application**

The **white** visitor's gown, **size XXXL**, is made of spunbonded nonwoven. The long gown is open in the back. It has long sleeves and an elastic band at the wrist, with ties at waist and collar. Due to the breathable material, the product is comfortable to wear. The white colour of the gown can be used for the assignment to specific sectors and for the identification "visitors". The use of the gown reduces the spread of germs through contaminated street clothes from the visitor to the patient in particularly sensitive areas (for example intensive care units).

### **Composition**

Polypropylene  
Sleeve cuff: polyester

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulatory Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

### **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### **Storage and Transport**

Dry and dustfree.

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.