



## BESUCHERKITTEL (Visitor's Gown)

REF 665705

### Product Description, Intended use, Application

The **green** visitor's gown, **size M**, is made of spunbonded nonwoven. The long gown is open in the back. It has long sleeves and hydrophilic knitted cuffs at the wrist, with ties at waist and collar. Due to the breathable material, the product is comfortable to wear. The green colour of the gown can be used for the assignment to specific sectors and for the identification "visitors". The use of the gown reduces the spread of germs through contaminated street clothes from the visitor to the patient in particularly sensitive areas (for example intensive care units).

### Composition

Polypropylene  
Sleeve cuff: polyester

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.