



NOBAMED Paul Danz AG

SCHUTZAUFLAGE

REF 667175

Product Description, Intended use, Application

The protection sheet, **size 80 x 175 cm**, is provided with a transparent plastic coating which is impermeable to liquids. The absorbent protection sheet rapidly absorbs liquids, while the foil coating on its underside prevents that the sheet is soaked. The product is reinforced with **20 blue threads** to increase its tensile strength and stability. The protection sheet is used in medical sectors and in nursing to protect stretchers, couches, beds, or tables from soiling, thus reducing the risk of contamination and counteracting the spread of germs. The protection sheet is for single use.

Composition

Cellulose, polyethylene, polyester

Contra-Indications

No contra-indications known.

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN 1041.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.