



**NOBAMED Paul Danz AG**

# SCHUTZAUFLAGE

**REF 667210**

## **Product Description, Intended use, Application**

The protection sheet, **size 80 x 210 cm**, is provided with a transparent plastic coating which is impermeable to liquids. The absorbent protection sheet rapidly absorbs liquids, while the foil coating on its underside prevents that the sheet is soaked. The product is reinforced with **20 blue threads** to increase its tensile strength and stability. The protection sheet is used in medical sectors and in nursing to protect stretchers, couches, beds, or tables from soiling, thus reducing the risk of contamination and counteracting the spread of germs. The protection sheet is for single use.

## **Composition**

Cellulose, polyethylene, polyester

## **Contra-Indications**

No contra-indications known.

The product should not be used in the case of a known allergy against the material.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Normative and Regulative Requirements, Common Standards**

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

## **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN 1041.

## **Storage and Transport**

To be stored in a dry and dust-free environment

## **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

## **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.