



NOBAMED Paul Danz AG

NOBA[®] OP-MASKE

REF 672032

Product Description and Purpose

The blue, 3-ply medical face mask with elastic ear loops has a filtration efficiency of $\geq 98\%$. Its differential pressure is very low. The microbiological filtration performance of all masks is very high, thus providing a safe protection for the patients and at the same time an additional protection for the user from splashes of contaminated liquid. The masks are equipped with a worked-in nose bow which can be matched to the dorsum of the nose. Edges and inner non-woven of the masks are extra soft. The masks are free from glass fiber and latex.

Composition

Polypropylene, latex-free elastomer thread, nose bow: stainless steel

Normative and Legal Requirements

NOBA[®] OP-Maske is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1. It meets the requirements of DIN EN 14683, type II R ("Medical face masks, Requirements and Test Methods").

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1-labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

