



**NOBAMED Paul Danz AG**

## **NOBA<sup>®</sup> OP-MASKE**

**REF 672047**

### **Product Description and Purpose**

The green, 3-ply medical face mask with tie-on bands and a **transparent visor** has a filtration efficiency of  $\geq 98\%$ . Its differential pressure is very low. The microbiological filtration performance of all masks is very high, thus providing a safe protection for the patients and at the same time an additional protection for the user from splashes of contaminated liquid. The masks are equipped with a worked-in nose bow which can be matched to the dorsum of the nose. Edges and inner non-woven of the masks are extra soft. The masks are free from glass fiber and latex.

### **Composition**

Polypropylene, Polyester  
Nose bow: stainless steel

### **Normative and Legal Requirements**

NOBA<sup>®</sup> OP-Maske is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1. It meets the requirements of DIN EN 14683, type II R ("Medical face masks, Requirements and Test Methods").

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1-labels on all its packaging.

### **Packaging**

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment

**The product bears following labels and symbols:**

