



NOBAMED Paul Danz AG

NOBA[®]-OP-Haube **Surgeon cap astronaut**

REF 673190

Product Description, Intended use, Application

The surgical caps (size M) are made of a spunlace nonwoven. The material is permeable to air and therefore breathable. By this an air This prevents heat accumulation, which increases the wearing comfort.

The surgical cap is used in the OR area to cover hair in a functional and safe way. It is provided with a hydrophilic sweatband in the front area.

Composition

Viscose, polyester, latex-free elastomer thread

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: Carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree, protect from sunlight.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.