



**NOBAMED Paul Danz AG**

## **NOBA<sup>®</sup>-OP-Haube**

### **Surgeon cap astronaut**

**REF 673200**

#### **Product Description, Intended use, Application**

The surgical caps (size L) are made of a spun-lace nonwoven. The material is permeable to air and therefore breathable. By this an air This prevents heat accumulation, which increases the wearing comfort.

The surgical cap is used in the OR area to cover hair in a functional and safe way. It is provided with a hydrophilic sweatband in the front area.

#### **Composition**

Viscose, polyester, latex-free elastomer thread

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### **Normative and Regulative Requirements, Common Standards**

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

#### **Packaging**

Primary packaging: folding box made of cellulose

Secondary packaging: Carton made of cellulose

#### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

#### **Storage and Transport**

Dry and dustfree, protect from sunlight.

#### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.