



NOBAMED Paul Danz AG

HAUBE

Product Description and Purpose

The caps are made of a spinning non-woven which is air-permeable and breathable. In this way a heat accumulation is prevented and the wearing comfort is increased. The caps are available in several different versions.

Composition

Polypropylene, elastomer thread (latex-free)

Normative and Legal Requirements

The cap is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

