



**NOBAMED Paul Danz AG**

## **DRAINAGE SET**

### **Product Description and Purpose**

The sterile, combined set is suited for the care of different drainages. It contains the following articles:

2 x non-woven slit swabs NOBATOP®S,  
5 cm x 5 cm  
2 x non-woven swabs  
NOBATOP® 8.5 cm x 5 cm  
1 x adhesive non-woven RUDAVLIES®  
10 cm x 10 cm

### **Composition**

Viscose, polyester, polyacrylate  
adhesive

### **Normative and Legal Requirements**

The drainage set is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH. It has DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: film-film  
packaging

Secondary packaging: folding box  
made of  
cellulose

Tertiary packaging: carton made of  
cellulose

### **Storage**

To be stored in a dry and dust-free environment, protected against direct sunlight.

**The product is bearing the CE marking, in compliance with the directive's provisions.**

