



NOBAMED®

Product Description and Purpose

The sterile pouch contains separately 2 + 3 plum size non-woven gauze balls. Suited for the disinfecting preparation and aftertreatment of angioaccesses.

Composition

Viscose, polyester, silicone

Normative and Legal Requirements

NOBAMED® is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper-film packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment.

The product is bearing the CE marking, in compliance with the directive's provisions

