



NOBAMED Paul Danz AG

WATTESTÄBCHEN

Cotton tipped on both sides

REF 974202

Product Description intended use, application

The cotton applicator, which is cotton tipped on both sides, is used as applicator for liquids and ointments, for sanitation and for collection of smear samples. For use on open wounds, the product must be sterilized with a validated sterilization method before (moist heat according to DIN EN ISO 17665 or ethylene oxide according to DIN EN ISO 11135).

Composition

Cotton, polypropylene

Note

Sterilize before use on open wounds.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

The cotton tipped applicator is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: Plastic box

Secondary packaging: Carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.